


Easy Setup Software




➔

Previous vs Current

Current workflow has all automated steps set up

Single Unified Control Interface



Easy Setup Software

EPIC Controller Dashboard

● READY ● BLOCKED ● UNKNOWN

Search controllers

combined_state_controller

Run

exo_controller_daniel

Run

exo_controller_daniel_debug

Run

hike_controller

Run

hip_controller_COMBINED

Run

hip_controller_COMBINED_flexible_devices

Run

hip_controller_COMBINED_noinsole

Run

hip_controller_COMBINED_noinsole_BO

Run

hip_controller_COMBINED_noinsole_BO_v2

Run

hip_controller_COMBINED_noinsole_BO_v3

Run

hip_controller_COMBINED_noinsole_BO_v4

Run

hip_controller_COMBINED_noinsole_BO_v5

Run

hip_controller_COMBINED_noinsole_BO_v6

Run

hip_controller_COMBINED_noinsole_BO_v7

Run

hip_controller_KEY

Run

hip_controller_noback-xs

Run

hip_controller_noback

Run

model_comparison_controller

Run

prueba

Run

SENSORS

Device ID	Type	Status
154134	IMU	connected
154135	IMU	connected
154136	IMU	connected
154137	IMU	connected
154138	IMU	connected

ACTUATORS

Motor ID	Type
2	AK80-9
4	AK80-9

SUMMARY

Jetson Hostsully: 192.168.50.220

XsensensorsOFF

Wi-FiCAREN_5G

Last updated2025-12-01 20:10:57

Flexible Controller - Devices

IMU IDs

Pelvis

✓ Select IMU

154134

154135

154136

154137

154138

Right Thigh

Select IMU

Right Shank

Select IMU

Left Shank

Select IMU

Motors

Left Hip

Motor ID

Select Motor

Type

—

Right Hip

Motor ID

Select Motor

Type

—

Reset Motors

☒ Use Xsensors (insoles)

Cancel

Run

Next Steps

- Explore on-device electronics packaging to eliminate need for separate backpack subassembly, beginning with integrating Jetson and PCB directly onto exoskeleton frame
- Streamline Easy Setup Software with current GUI and package into easily accessible web app
- Continue testing don/doff prototype and build final assembly
- Manufacture and test new backplates